Titus

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Contents

	Page
For Your Safety and that of Your Patients	3
Recommendations on handling oxygen and compressed-gas cylinders.....	4
Intended Use	5
Monitoring	6
Description of Models	7
Assembling Titus	10
Gas supply.....	10
Assembling the circle system.....	13
Measuring the expiratory minute volume.....	13
Connecting hoses.....	14
Compressed gas supply for Ventilog.....	15
Combination with Ventilog 2.....	16
Combination with Ventilog 3.....	17
Fitting Vapor.....	18
Monitoring.....	19
Connecting to anaesthetic gas scavenging system.....	22
Manual ventilation bag.....	23
Equipotential bonding.....	23
For emergency use only.....	24
Titus MRI for Use with MRI	25
Equipment Check	31
Checklist	34
Operation	37
Fresh gas supply.....	37
Operating modes.....	39
For emergency use only.....	41
Checking use.....	42
What to do if gas supply fails.....	42
Additional instructions.....	43
End of use.....	44
Care	45
Stripping down machine.....	45
Magill inhalation device.....	45
Care list for Titus anaesthetic apparatus.....	48
Cleaning and disinfection.....	50
Sterilisation.....	51
Service/Maintenance intervals.....	51
Testing Readiness for Operation	52
Appendix	58
Functional Diagrams.....	58
Safety Features.....	61
Technical Data.....	63

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For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observance of these instructions. The apparatus is only to be used for purposes specified here.

Maintenance

The apparatus must be inspected¹⁾ and serviced¹⁾ by experts at regular 6 month intervals (and a record kept). Repairs¹⁾ and general overhaul of the device may only be carried out by experts. We recommend that you take out a service contract with DrägerService and ensure that all repairs are carried out by them. Only original Dräger parts must be used for maintenance¹⁾. Observe the chapter on "Maintenance Intervals".

Accessories

Only use the equipment itemized on the list of accessories.

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1)
 Inspection = examination of actual condition
 Service = measures to maintain specified condition
 Repair = measures to restore specified condition
 Maintenance = inspection, service and, if applicable, repair

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Dräger Medical AG & Co. KGaA

Safe connection with other electrical equipment

Electrical connections to equipment which is not listed in these book should only be made following consultations with the respective manufacturers or an expert.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorised by DrägerService or if the apparatus is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendation given above.

The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

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Recommendations on handling oxygen and compressed-gas cylinders

- Oxygen cylinder valves and oxygen pressure reducers must not be oiled or greased, or touched with oily fingers.

Danger of explosion!

- Do not store oxygen cylinders with substances which ignite easily.
- Protect filled cylinders from direct sources of heat (do not store in direct sunlight or close to a radiator).
- Valve handwheels should only be operated manually. Do not use tools. Cylinder valves are precision components which can easily be damaged when force is applied.
- Leaking or stiff cylinder valves should be repaired by trained personnel.

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Intended Use

Inhalation anaesthetic machine with continuous fresh gas flow. Fresh gas flow via a flowmeter block for O₂ and N₂O (optional: also for air) and via an anaesthetic vaporiser.

The machine is intended for mobile use for First Aid at the scene of an accident^{*)} and during transport (gas supply from a portable cylinder pack) as well as for use in medical environments, such as surgeries and hospitals; in these situations, gas is normally supplied from a pipeline system (PS). Supply is possible from 3-litre or 11-litre cylinders (O₂ and N₂O) if Titus is placed on a specially equipped trolley.

In a specific version, Titus can be used in magnetic resonance imaging (see page 25 »Titus MRI«).

Titus can be fitted with the ISO circle system 9 or a semi-open anaesthetic system or alternatively a Magill inhalation device.

Optional:

ISO Infant Circle System for small children.

Instructions for use

Only not inflammable anaesthetic agents may be used according to EN 740! Danger of fire!

As this machine is not licensed for the use of inflammable anaesthetic agents (ether, cyclopropane or similar agents), the use of antistatic (conductive) respiratory tubes and face masks is not necessary.

Conductive respiratory tubes and face masks can cause burns at high-frequency surgery and in this case are not recommended for this machine.

Additional electrical appliances that are latched onto the top, have to be connected to the main machine by means of a potential compensation line.

Use the machine only under supervision of qualified medical personnel in order to react immediately and appropriately in case of malfunction.

^{*)} At the scene of an accident, only machines with a fixed anaesthetic vaporiser may be used (see Fixed Vapor p. 9)

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Monitoring

For anaesthetic workstations, the European Standard (EN 740)* requires monitoring and alarm in order to protect against dangerous initial conditions:

- airway pressure measuring with alarm upper limit (alarm lower limit is recommended).
- expiratory tidal volume measuring (alarm lower limit is recommended).
- inspiratory oxygen concentration with alarm lower and upper limit.
- measuring of the concentration of anaesthetic agents with alarm lower and upper limit.
- CO₂ concentration measuring with alarm lower and upper limit.
- checking of disconnection and stenosis in case of automatic ventilation.

These facilities are not part of the Titus. Therefore adequate additional appliances have to be provided. Relevant interfaces are described in the technical data. Is the Titus deployed as a First Aid unit, the requirements above are only recommendations.

For reason of safety, Titus with integrated auxiliary compressed-air unit must be used with an inspiratory O₂ monitor.

If monitoring equipment is used, the user must check that the equipment conforms to legal requirements, and reliably monitors the efficiency of the anaesthetic ventilator.

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*) Anaesthetic workstations and their modules
– particular requirements –

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

Description of Models

Basic Units

Items of equipment in brackets () are not part of the basic equipment and are supplied only on special order.

- Portable anaesthetic machine to be hooked on a wall rail on the patient's bed or on a (trolley).
- Connections for supplying compressed gases O_2 , N_2O (and Air) from a pipeline system – PS – (or from gas cylinders).
- Flowmeter block for measuring flow of O_2 and N_2O (and Air).

On special order

- Trolley without holders or with holders for two 3 or 11-litre O_2 and N_2O cylinders. Cylinders, pressure reducers and connecting hoses must be ordered separately. Monitors may be mounted on the upper shelf of the trolley. The ventilator can be mounted with an additional baseplate.
- Titus-Air with additional probe for compressed air (see top picture).
- Only O_2 and N_2O can be supplied from cylinders – not Air. Different connecting hoses are required for a supply from cylinders than for a supply from a pipeline system (PS).
- Air and N_2O can only be dosed alternately via one and the same flowmeter tube (see also the diagram on page 60). N_2O or Air is selected manually using the switch below the flow valve.

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The following basic models are available.

Titus Low-Flow	Vapor-Stecksystem	mit S-ORC	M 33 041	Plug-in Vapor	with S-ORC
Titus Standard	Vapor-Stecksystem	mit S-ORC	M 33 042	Plug-in Vapor	with S-ORC
Titus Standard	Vapor-Festanschluß	mit S-ORC	M 33 043	Fixed Vapor	with S-ORC
Titus Standard*	Vapor-Stecksystem	ohne S-ORC*	M 33 044	Plug-in Vapor	without S-ORC*
Titus Standard*	Vapor-Festanschluß	ohne S-ORC*	M 33 045	Fixed Vapor	without S-ORC*
Titus NMR	Vapor-Stecksystem	mit S-ORC	M 33 046	Plug-in Vapor	with S-ORC
Titus Asia**	Vapor-Stecksystem	mit S-ORC	M 32 947	Plug-in Vapor	with S-ORC
Titus Asia**	Vapor-Festanschluß	mit S-ORC	M 32 948	Fixed Vapor	with S-ORC
Titus NMR Asia**	Vapor-Stecksystem	mit S-ORC	M 33 058	Plug-in Vapor	with S-ORC

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Machines without ORC may only be used in countries where ORC is not a legal requirement. These machines are not designed with the CE-Label.

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Description of safety features:

- O₂ shortage signal (in all basic units)
- S-ORC and
- N₂O cut-off

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*) In the area of validity of the European Standard not licensed.

**) Contrary to the European Standard, these machines have national colour codes and flowmeter tubes, which are positioned the other way round. No sale within the European economic area.

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
**Brief description****1 Plug-in Vapor:**

The Dräger-Vapor® * anaesthetic vaporiser can be unlocked with a lever and detached from the anaesthetic machine. This enables a quick change-over to a different anaesthetic agent.

Titus with plug-in Vapor must not be used for mobile applications (such as First Aid at emergency)!

2 Fixed Vapor:

The Vapor is fixed to the anaesthetic machine, and can only be changed by DrägerService.

Titus with a fixed Vapor is suitable for mobile use at emergency.

3 Sensitive Oxygen Ratio Controller S-ORC:

Flow-dependent control to maintain a minimum O₂ concentration (21 vol.%) in the fresh gas.

Machines with S-ORC are marked with an »S-ORC« label on the machine.

*) In this case it is either a Vapor 19.n or a Vapor 2000, in following only called as Vapor.

Attention: Using a Vapor 2000, the Titus has to be modified accordingly.

Modification: The rubber nozzle behind the Vapor must be replaced with a thin one. So that the Vapor can stand vertically.

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- Beide Dosierventile schließen.

Assembling Titus

The machine must not be tilted excessively!


The Vapor 19.n may be tilted as far as 45° without any deviation in concentration (until 30° at Vapor 2000).

During transport, when the Vapor may be tilted by more than 45°, it must be emptied of anaesthetic agent (only for Vapor 19.n). Adhere strictly to the instructions for use of the Vapor.

If the anaesthetic machine is placed on a sloping surface, the flowmeter may show greater variation than is specified under »Technical Data«.

To assemble, either

- hook the machine onto a wall rail or
- place it on a level, horizontal base (table or trolley)
- or
- hook it over the oval-shaped rail on the special Titus trolley.

Gas supply

The recommendations for »Air« given in brackets () apply only if the anaesthetic machine has been supplied for air as a third gas.

This is only possible if Titus is fitted with the appropriate compressed gas connection.

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
Supply from a pipeline system (PS)

- Screw hoses onto the connectors for O_2 and N_2O .
- (Screw hose onto the air connector on Titus – male thread on Titus.)
- Connect the probes on the hoses to the wall terminal units of the pipeline system.

O_2 and N_2O supplies from cylinders

- Place the 11-litre cylinders complete with cylinder jackets or the 3-litre cylinders in the holders on the trolley and secure with leather straps.
- Screw pressure reducers to cylinder valves.
- Connect hoses between pressure reducers and gas inlets on Titus.

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
O₂ and N₂O supplies from portable cylinder pack

for supply during mobile use and as an emergency supply in hospitals.

The portable cylinder pack – about 10 kg – contains one 2-litre cylinder each of O₂ and N₂O.

The cylinder pack can be hooked onto a wall rail or onto a bed, or placed on the floor.

- Supply gas to anaesthetic machine – either

a) via hoses with probes for the pipeline system or

b) via hoses with threaded connectors for pressure reducers. In this case, unscrew the terminal units from the cylinder pack.

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Assembling the Circle System

- Prepare and assemble components following »ISO Circle System 9 Instructions for Use«.
- Fit circle system arm over pin on anaesthetic machine, and secure with knurled screw.

- 1 Fill CO₂ absorber** with fresh soda lime (Drägersorb® 800) and fit to the input cone on the circle system arm.
- 2 Fit inspiratory valve** to the input cone on the CO₂ absorber; turn nozzle to the left.


For emergency mobile use:

- Connect airway pressure gauge to the circle system.

Measurement of expiratory minute volume

With PM 8050 or alternatively with the PM 8060 or Minute Volumeter.

- 3 Connect Minute Volumeter** to the circle system. Alternatively:
- 4 Screw on the sensor housing with flow sensor** for patient monitor PM 8050/8060 (see Instructions for Use "Patient Monitor PM 8050/8060").
- 5 Connect expiratory valve** to Minute Volumeter. Alternatively: to sensor housing.

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
Connecting hoses

Fit one breathing hose each to

- 1 inspiratory valve and
- 2 expiratory valve and connect to
- 3 Y-piece.
- Fit
- 4 mask or catheter connector to Y-piece, depending on the patient's needs.

Optional:

- Use microbial filter on inspiratory or expiratory valve.

Follow »Microbial Filter Instructions for Use«.

- 5 Connect breathing bag to
- 6 connector and
- 7 breathing hose when no automatic ventilator is in use. Otherwise: see page 16.
- 8 Fit fresh gas hose to anaesthetic machine and to circle system arm.

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
Compressed gas supply for Ventilog

Place ventilator beside Titus, preferably on a bracket mounted on the same wall rail as Titus.

Other options:

- Titus on trolley:
on separate latching plate of the trolley.
- Separate trolley with latching plate.
- Work surface.

Ventilog **must not be placed above** electrically-operated machines – for instance on a monitor. Danger of fire from the emitted pure oxygen or oxygen-enriched air. The reverse arrangement is permitted.

With a supply from the pipeline system use a medical air distributor or oxygen distributor with suitable connecting hose.


- PS-supply – Air as a third gas is available – Titus is additionally fitted with an air-connector.
- PS-supply – only O₂ is available – Titus is not fitted with an air connector (or air-connector remains open).

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Combination with Ventilog 2

Connecting Ventilog 2 with Circle System 9

Follow the Ventilog 2 Instructions for Use!


- 1 Screw the pneumatic changeover valve to circle system arm –
- 2 Connect the 1.1 m breathing hose to the Ventilog 2 anaesthetic ventilator and connect to the changeover valve –
- 3 Connect the control line – one nozzle to the front panel of the Ventilog 2 anaesthetic ventilator and the other to the rear of the changeover valve –
- 4 Connect the fresh gas hose –
- 5 Connect the breathing bag with the breathing hose to the changeover valve – hang up the breathing bag.

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Combination with Ventillog 3

Connecting Ventillog 3 with Circle System 9

Follow the Ventillog 3 Instructions for Use!


- 1 Screw the FGE pneumatic changeover valve to the circle system support arm –
- 2 Connect the 1.1 m breathing hose to the Ventillog 3 anaesthetic ventilator and to the FGE changeover valve –
- 3 Connect the fresh gas hose – from Titus to the FGE changeover valve –
- 4 Connect the breathing bag with the breathing hose to the FGE changeover valve – hook up the breathing bag –
- 5 Connect the control line – to the back of the Ventillog 3 anaesthetic ventilator and to the FGE changeover valve –
- 6 Connect the control line for PEEP – back of Ventillog 3 – and to the FGE changeover valve – press the plugs into the sockets until they click into place.

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Fitting Vapor

(only for Titus with plug-in Vapor)

- Follow »Vapor 19.n or Vapor 2000 Instructions for Use«!
- 1** The mounting on the anaesthetic machine holds the Vapor.
- Make sure that both O-rings are in place and undamaged.
- 2** Swing the locking lever on the Vapor to the front.
- 3** Fit the Vapor plug-in adaptor over the two studs on the
- mounting – must be level.
- Swing the locking lever to the left – as far as it will go – the Vapor is locked.

If no Vapor is fitted, the holes in the studs will be sealed, and there will be an automatic connection between flowmeter block and common gas outlet.

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
Monitoring

Measurement of expiratory minute volume, airway pressure, inspiratory O₂ concentration, anaesthetic concentration in the fresh gas and the inspiratory and expiratory CO₂ concentration – e.g. with monitors PM 8050 or PM 8060.

In the case of automatic ventilation, an alarm system must be provided for disconnection and stenosis.

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Mounting the monitor

- 1 Tilt the monitor forwards and insert the front latch in the cover plate slot –
- 2 Insert the rear latch in the slot and tighten with the knurled screw on the rear.

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
Installing and connecting PM 8050/PM 8060 sensors

- 1 T-piece: screw the test line tightly to the filter and water trap –
- 2 Optional: temperature sensor in the inspiratory side of the Y-piece –
- 3 Fit the O₂ sensor for inspiratory O₂ measurement – (not required for the PM 8060)
- 4 Fit the pressure measuring line for inspiratory airway pressure measurement –
- 5 Flow sensor for expiratory volume measurement – electrical connection –

Connect the sensor lines to the back panel of the monitor in accordance with the "PM 8050" or "PM 8060" Instructions for Use.

- Connect the mains cable, firstly to the breathing gas monitor and then in the electrical socket.

Monitoring the inspiratory O₂ concentration

e.g. with Oxydig

- 6 Place the measuring device on the holder. No mains connection required. Device operated by batteries.
- 7 Connect the O₂ sensor with Oxydig.
- 8 Fit the O₂ sensor.

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Monitoring the expiratory and inspiratory CO₂ concentration respectively anaesthetic gas

e.g. with PM 8050

alternatively:

PM 8060.

- Hang Titus on the wall rail.
- Place the PM 8050 monitor on a shelf right to the Titus or latch it onto the base plate of the Titus trolley. (If necessary, remove both baseplate in order to latch).

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Connecting sampling hose

For sampling to measure CO₂, anaesthetic gas and O₂:

- 1 Screw the filter into the T-piece.
- 2 Fit T-piece of the patient's connector into the Y-piece - filter at the top to avoid occlusion because of moisture.

3. Screw sampling line firmly onto the filter and water trap.


Use only the original sampling line; other lines may distort the machine's technical data.

4. A Y-piece with Luer Lock can be used instead of a T-piece and the filter.

Do not allow alcohol or alcohol-based agents to get into the sampling line!

Alcohol would lead to incorrect concentration measurements.

- Follow the Instructions for Use of the monitor.

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
Connecting to anaesthetic gas scavenging system (AGSS)

- Fix the AGSS with the relevant kit.
- Connect transfer hoses to the over pressure valve and (if available) Ventillog.
- Connect the scavenging hose to the scavenging system.

Anaesthetic gas scavenging

According to the European Standard EN 740, an Anaesthetic Workstation can be used only with the AGSS. The following systems below do not fulfill the requirements.

Connecting to anaesthetic gas scavenging

- Fit scavenging hose to the nozzle below the relief valve (picture: view from the back).
- Connect the probe on the scavenging hose to the terminal unit on the scavenging system only when anaesthesia is about to commence.

If Ventillog is being used:

- also connect hose to the nozzle on the back of Ventillog.

If no scavenging system is available:

- Fit an »anaesthetic filter set 2« to each unit (corrugated hose fitted with Anaesthetic Filter 633).

Attention:


The filters saturate and have to be changed in time. There is no N_2O -binding by the filters. The therefore developed dangerous influence can be reduced buy increased ventilation of the room as mentioned above.

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Manual ventilation Device (ventilation bag)

such as Resutator 2000 for emergency ventilation.

Place the ventilation device close to Titus for example on a wall rail or trolley.


If a life-threatening malfunction is detected at the ventilator, the patient must be ventilated with a separate manual ventilation device or a breathing bag.

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Equipotential bonding

for intracranial and intracardial operations.

If the anaesthetic machine is mounted on a trolley:


- Establish the earth bonding with a suitable cable between the bolt on the back of the trolley and the equipotential bonding contact of the room.
- Also install equipotential bonding for the breathing gas monitors!

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For emergency use only

Connection of a Magill Inhalation Device.

The picture shows the Magill Inhalation Device ISO- M 26130.

- 1** Screw the Magill connector onto the ISO fresh-gas outlet.
- 2** Attach the 2.3 L reservoir bag.
- 3** Connect on a 1.1 m breathing hose.
- 4** Push the mask nozzle into the breathing hose.
- 5** Push the mask onto the mask nozzle.
- 6** Turn the lever to the vertical position.

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Checking the machine

immediately before use.

A checklist must be fixed to the anaesthetic machine.

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for use with MRI

Identification

Two components of the anaesthetic system carry the above label, which refers to the modification of these components for use in MRI-scanning:

- Titus basic unit (on the front)
- Ventillog 2 or Ventillog 3 (on the back).

The Titus trolley is suitable for the use in a MRI suite and this label is not attached on it.

These components can also be used for other operating theatre applications.

Approved accessories

The Circle System 9 (optionally the infant circle system set) with extended ventilation hoses (3 ... 4 m) has been approved as ventilation system.

The following have been approved as monitors for use in MRI:

- Volumeter 3000 and 2000 K (for children),
- Oxydig (O₂ monitoring)
- PM 8050 NMR (respiratory gas monitor)

It is not approved to operate the Titus MRI with gas cylinders in a MRI suite.

The examination room must be equipped with a pipeline system, an anaesthetic scavenging system and a 230 V AC main power supply.

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Safety regulations and recommendations

for use in MRI scanning

- Compressed gas cylinders are not allowed in examination rooms.
- Always lock castors on trolley.
- Only accessories approved for use in scanning may be used.
- Batteries must not be changed on monitors (such as Oxydig) in the examination room.
- Vapor must always be changed away from the MRI examination room.
- Do not attach an Oxydig O₂ measuring device to the Titus housing if the Titus is placed in a field higher than 30 mTesla.
- Extended breathing hoses (3 ... 4 m) result in increased compliance of the breathing system.
When using ventilation hoses for adults, increasing the set volume (V_T) by about 5 ... 10 % may be appropriate within a range of 100 ... 600 mL tidal volume (V_T).
- When ventilating children: do not use ventilating hoses with metal inserts. Metal inserts result in a distortion of the MRI picture.
- Silicone hoses used in anaesthesia require special cleansing treatment cycles. Follow the manufacturer's instructions.
- Do not use any intubation accessories made of ferromagnetic metal in the MRI examination room, such as endotracheal catheters, mouth tubes, laryngoscopes etc..
- A Half-yearly functional inspection at a MRI scanner is required. The anaesthetic machine and the monitoring equipment, including sensor extension cables, must also be checked.

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
Installation of anaesthetic machine and monitoring in the MRI examination room

The Titus MRI inhalation anaesthetic machine and the monitoring accessories listed previously have been tested in combination with MRI scanners up to 1.5 Tesla.

Field intensities vary according to the type of scanning unit, and so extra care and attention are required during installation and operation of an anaesthetic machine with accessories.

Installing anaesthetic machine

The modified Titus MRI anaesthetic machine must not be operated in a magnetic flux density greater than 70 mTesla.

Installing the Titus MRI with trolley and fixed PM 8050 MRI

PM 8050 is firmly installed onto the Titus NMR with trolley.

- Place the trolley of the Titus MRI with the PM 8050 right or left of the tube.
- Keep the machine outside the 20mT-line.
- **Position the machine outside the 20mT-line, however not closer than 1.5 m from the center of the tube opening.**

Use antimagnetic tape measure.

- Lock castors of the trolley.
- Mark the position of the trolley, e.g. by a yellow-black tape on the floor and put it down on the installation plan of the machine.

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Installing monitors

All monitors – except the Volumeter 3000/2000 K – must be operated at a greater distance from the scanner than an anaesthetic machine with circle system, as they are influenced by magnetic fields and high frequency (picture distortion, malfunctions).

The Oxydig measuring and warning device can be installed and operated outside the 30 mTesla area. It is, however, recommended that this device is placed next to the other monitors, outside the 20 mTesla area.

The hoses and cables to the anaesthetic machine must be laid along the floor and secured with a cable cover which will withstand squashing, kinking and being tripped over.

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Documentation of machine hand-over

- Before the machine is handed over, DrägerService carry out a function check on-site.
- The start-up* and hand-over on site are documented and signed by a Dräger representative and the user of the unit. The record contains details of the distances from the scanner, in accordance with the scanner manufacturer's installation plan.

This document is an integral part of the Titus MRI machine approval and the original, or a copy, must be stored in the machine log.

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
Betreiber:

Standorte

Hand-over record and machine installation plan

Titus MRI with PM 8050 MRI for MRI scanning

Manufacturer of anaesthetic machine: Dräger Medical AG & Co. KGaA, Lübeck

Machine type: Titus MRI

Manufacturer of MRI scanner: _____

Machine type: _____

Place of installation: _____


Date of installation of anaesthetic machine: _____

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Equipment Check

A plastic checklist is supplied with the machine.

- By deleting and adding, adapt the checklist in the Instructions for Use (pages 34/35) to cover the types and variants of the equipment and accessories used.

This checklist will then contain all the checks which have to be carried out on the equipment before each use.

- In addition, add the following:
Model designation and serial number of the machine.
- All permanent entries on the plastic checklist must be made in waterproof »permanent« felt-tip pen.
- Attach the checklist to the machine with the bead chain.

Carry out equipment check immediately before each use!

- Enter all checks on the plastic checklist -
– in pencil -
– in the ACTUAL column -
– including date -
– and sign.
- Rub out all entries before the next check.
- Do not wipe the plastic checklist with cleaning agents, disinfectants, alcohol or similar solvents, because the waterproof felt-tip pen may be removed by these agents.

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Name:
Datum:

Check List for Titus, Sulla 909, Trajan 808, Livius M

All checks must be carried out before the equipment is used.

The person who carries out the checks must be fully conversant with the Instructions for Use .

Checks marked with a P must be carried out before use on each patient.

Anaesthesia machine,
Serial No:

Ventilator, Serial No:

Monitor, Serial No:

Tick each function when checks have been satisfactorily completed.

Pre-conditions

- Inspection intervals for machines and accessories not exceeded, eg. see check labels ☐
- Machine fully assembled ☐
- 2.3 L breathing bag fitted to patient connector as test lung on the Y-piece or filter, whichever is being used ☐
- Anaesthetic gas scavenging hose connected to terminal unit and to waste gas sockets on machine ☐
- Monitor (O₂, P, V, CO₂, anaesthetic agent) switched on, self-test carried out satisfactorily, sensor calibrated, sampling line for gas monitoring connected to Luer lock on the Y-piece, correct anaesthetic agent selected ☐
- Devapor (when being used) switched on ☐

Tests

Gas supply

- Open reserve gas cylinders (when present):
O₂ pressure more than 50bar, ☐
N₂O more than 30bar. ☐
- Close reserve gas cylinders.

Medical gas pipeline (when present):

Probes for O₂, AIR, N₂O connected to terminal units ☐

Flowmeters

for Sulla 909, Trajan 808 and Livius M:

- Turn AIR/N₂O selector switch to N₂O,
- Open O₂ flow-control valve fully:
O₂ flow present. ☐
- Open N₂O flow-control valve fully :
N₂O flow present. ☐
- Turn AIR/N₂O selector switch to AIR:
No N₂O flow present. ☐
- Open AIR flow-control valve fully:
AIR flow present. ☐
- Close AIR flow-control valve,
- Turn AIR/N₂O selector switch to N₂O.

for Titus:

- Turn AIR/N₂O selector switch to N₂O,
- Open O₂ flow-control valve fully:
O₂ flow present. ☐
- Open N₂O/AIR flow-control valve fully:
Flow present in N₂O/AIR flowmeter. ☐
- Turn AIR/N₂O selector switch to AIR:
Flow present in N₂O/AIR flowmeter. ☐

- Disconnect Air probe:

No AIR flow present in N₂O/AIR flowmeter. ☐

- Connect AIR probe,
- Turn AIR/N₂O selector switch to N₂O.

S-ORC Function

Shut-down

- O₂ flow-control valve:

No N₂O flow. ☐

Control

- Set O₂ flow -control valve to 1.5 L/min.:

N₂O flow 3 to 5 L/min. ☐

- Close O₂ and N₂O flow-control valves.

Vapor 19.n / Vapor 2000

Fastening

- P Latched down firmly and set vertically ☐

Handwheel

- In zero position and engaged ☐
- P Filling level between min. and max. ☐

Interlock (when present)

- P Locking function OK ☐

Key-indexed filling system

Sealing key or pin inserted and close tightly

or

Quik Fil system or Funnel filling system

Locking screw tight ☐

Devapor (when present)

Fastening

- P Latched down firmly and resting evenly on surface ☐

Handwheel

- In zero position and engaged ☐
- P Filling level between half full and max. ☐
- "Operation" light lit ☐

Key-indexed filling system

Filling opening locked shut and pointing downwards ☐

Circle system

- Valve disks, valve supports in proper condition ☐
- Retaining nut screwed tight ☐

Soda lime

- Replaced or no discolouration ☐

Leak-tightness

- P** ● Switch-pressure limiting valve to "MAN" and set to 70 mbar,
- For Sulla 909, Trajan 808 und Livius M:**
- P** ● Set O₂ flow-control valve to 0.12 L/min.,
- P** ● Press O₂ flush; let pressure build up to ca. 30 mbar:
- P** ● Pressure remains constant or rises. ☐
- Switch pressure-limiting valve to "SPONT": Pressure drops.
 - Switch pressure-limiting valve to "MAN",
 - Set Vapor to lowest concentration shown on handwheel,
 - Press O₂ flush; let pressure build up to ca. 30 mbar: Pressure remains constant or rises. ☐
 - Reset Vapor to 0.
 - Release pressure: Push lever on pressure-limiting valve.
- For Titus:**
- P** ● Remove test lung,
- P** ● Detach hose from breathing bag and connect to patient connector or to filter.
- P** ● Press O₂ flush; let pressure build up to ca. 30 mbar:
- P** ● Pressure change from 30 mbar to 20 mbar takes at least 16 seconds. ☐
- Switch pressure-limiting valve to "SPONT": Pressure drops. ☐
 - Switch pressure-limiting valve to "MAN",
 - Set Vapor to lowest concentration shown on handwheel,
 - Press O₂ flush; let pressure build up to ca. 30 mbar: Pressure change from 30 mbar to 20 mbar takes at least 16 seconds. ☐
- P** ● Detach hose from Y-piece or filter and fit to breathing bag,
- P** ● Keep Y-piece sealed,
- P** ● Press O₂ flush; let pressure build up to ca. 30 mbar: Breathing bag does not empty when squeezed. ☐
- Reset Vapor to 0.
- P** ● Reconnect test lung.

Inspiratory and expiratory valves

- Set pressure-limiting valve to 30 mbar,
- Press O₂ flush; let pressure build up to ca. 30 mbar: Breathing bag for manual ventilation fills. ☐
- O₂ concentration on monitor more than 80 Vol.%. ☐
- Valve disc in insp. valve moves when breathing bag is squeezed during manual ventilation. ☐
- Valve disc in exp. valve moves when breathing bag is released. ☐
- Test lung on patient connection fills and empties rhythmically during manual ventilation. ☐
- No resistance during manual ventilation which might indicate stenosis. ☐

Pressure-limiting valve

- Set O₂ flow-control valve to 5 L/min.: Display on monitor: 25 to 35 mbar. ☐

Ventilator

- Set stroke volume to 600 mL,
- Set frequency to 10/min.,
- Set PEEP to "12 o'clock position",
- Set other parameters to "green dot":

Display on monitor should be:

- Stroke volume 450 to 750 mL ☐
- Frequency 8 to 12/min ☐
- PEEP present. ☐

- Switch off ventilator,
- Set PEEP to 0,
- Close O₂ flow-control valve,
- Switch monitor to Standby.

Secretion suction

- P** ● Switch on suction,
- P** ● Kink suction hose: Display on pressure gauge between – 0.5 and – 0.9 bar. ☐
- P** ● Set negative pressure to suit patient's needs and switch off suction.

Breathing bag for emergency ventilation

- Present and ready for use ☐

Anaesthetic gas scavenging

- Float in AGS flowmeter between minimum and maximum ☐

Before connecting to patient

- P** ● Pre-set tidal volume and frequency to suit patient's needs.
- P** ● Set alarm limits for O₂, P, V, CO₂ and anaesthetic agent to suit patient's needs.
- P** ● Set pressure-limiting valve to suit patient's needs.
- P** ● Set monitor to anaesthetic agent being used.
- P** ● For monitors with sidestream gas measurement empty water trap.
- P** ● When there is electrical contact to patient: Establish potential equalisation between machines.

If any check cannot be carried out satisfactorily, the machine must not be used.

Machine check carried out

Name:

Date:

Was/What	Wie/How	Soll/Desired	Ist/Act.

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
Operation

Before operating the anaesthetic machine, make sure that:

- Readiness for operation has been checked (pages 52 to 57),
- Checks have been carried out according to the checklist (pages 31, 34, 35)
- Instructions for Use of all accessories, such as the Circle System, Vapor, anaesthetic machines, monitors, etc. are understood.

Fresh gas supply

Only for Titus, which has air supplied as a third gas:

- 1 Move selector switch to the desired position:
»N₂O« (O₂ + N₂O can be metered)
or
»Air« (O₂ + Air can be metered)
For Air/O₂ mixtures, see table on the following page.

Depending on the setting of the selector switch, N₂O or Air flows through the right-hand flowmeter tube.

If an additional compressed air device is integrated in the system, inspiratory O₂ measurement is required.

- 2 Set mixing ratio and flow of fresh gas with flow control valves.
Only when fresh gas flow has been set:
 - 3 Switch Vapor on:
keep white button pressed and set desired concentration with handwheel
or
switch Vapor off: turn handwheel to zero, with button for zero-point lock engaged.
- Connect Y-piece of the completely assembled circle system to patient via mask or tube.

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Variable O₂ concentrations in
Air/O₂ mixtures

Vol. % O ₂	Flow in L/min (Air + O ₂)													
	1		2		4		6		8		10		12	
	Air	O ₂	Air	O ₂	Air	O ₂	Air	O ₂	Air	O ₂	Air	O ₂	Air	O ₂
21	1.0	0	2.0	0	4.0	0	6.0	0	8.0	0	10.0	0	12.0	0
30	0.9	0.1	1.8	0.2	3.5	0.5	5.3	0.7	7.1	0.9	8.9	1.1	10.6	1.4
40	0.8	0.2	1.5	0.5	3.0	1.0	4.6	1.4	6.1	1.9	7.6	2.4	9.1	2.9
50	0.6	0.4	1.3	0.7	2.5	1.5	3.8	2.2	5.1	2.9	6.3	3.7	7.6	4.4
60	0.5	0.5	1.0	1.0	2.0	2.0	3.0	3.0	4.0	4.0	5.0	5.0	6.0	6.0
70	0.4	0.6	0.8	1.2	1.5	2.5	2.3	3.7	3.0	5.0	3.8	6.2	4.6	7.4
80	0.3	0.7	0.5	1.5	1.0	3.0	1.5	4.5	2.0	6.0	2.5	7.5	3.0	9.0
90	0.1	0.9	0.3	1.7	0.5	3.5	0.8	5.2	1.0	7.0	1.3	8.7		
100	0	1.0	0	2.0	0	4.0	0	6.0	0	8.0	0	10.0		

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Any flow deviations during use with »Air« can be found in the graph given in the Technical Data (page 65).

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
Spontaneous breathing

with Circle System 9

- 1 Turn lever of pressure limiting valve to the »SPONT« position; the patient can now breathe out freely via the check valve.
- 2 Fill the breathing bag sufficiently – supplies ready-metered fresh gas for spontaneous breathing.

Manual ventilation


with Circle System 9

- 3 Turn lever on the pressure limiting valve to the »MAN« position, and set the maximum airway pressure.
Ventilate manually with
- 4 breathing bag. Make sure it is full enough. Do not kink or squash hose.

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Automatic ventilation

with Ventillog 2 or Ventillog 3 and Circle System 9. Shown with Ventillog 2.

- 1 Turn the lever of the pressure limiting valve to »MAN« and set the pressure limit on the scale and
- 2 Set the »Ventilog« anaesthetic machine for automatic ventilation.

- the pneumatic switching valve switches to »Autom.«.
- Set the ventilation pattern on the Ventillog, according to the relevant Instructions for Use.

To monitor ventilation:


- Switch on patient monitor PM 8050. Set according to relevant Instructions for Use.

Changeover from automatic ventilation to spontaneous breathing or manual ventilation.

- Switch off Ventillog: turn rotary knob to »0« – the pneumatic switching valve switches automatically to »SPONT / MAN« –
- Continue as for »Spontaneous breathing« or »Manual ventilation«.

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Using O₂ flush (bypass)

- Press button – O₂ at about 55 L/min flows into the circle system to the patient, bypassing the flowmeter block and Vapor.
- Make sure that airway pressure does not increase unduly!

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Magill Inhalation System*

Before connecting to patient:

- 1 Turn the lever to the vertical position – rebreathing (exhalation) into the breathing bag is possible.
- 2 Set the O₂/N₂O mixing ratio and the flow of the fresh gas using the flow control valve. Fresh gas volume = minute volume. Smaller volumes of fresh gas result in partial rebreathing** and shortage of breathing gas, and the breathing bag is emptied (danger of suffocation!). Increase the fresh gas flow.
- 3 Switch on Vapor: Keep the white button pressed and set the concentration with the handwheel. For further information, see the »Vapor Instructions for Use«.
- 4 Turn the pressure relief valve anti-clockwise as far as it will go = this reduces any over pressure in the system to a minimum.

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If no rebreathing (or only partial rebreathing) is required:

- 1 Turn lever to a horizontal (or slanted) position.
- Increase volume of fresh gas as a result!

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* This device relies on the patient's own breathing (spontaneous breathing).

** Return breathing causes increased CO₂ concentration in the breathing gas.

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Checking use – of CO₂ absorbers and filters

CO₂ absorber

See »Circle System 9« Instructions for Use

Microbial filters, anaesthetic filters*

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What to do if gas supply fails

The following applies to all Titus machines – whether
fitted with S-ORC or not, or whether in a version with or
without air:

- If there is an O₂ shortage or failure:
acoustic alarm, for at least 7 seconds.
- If there is an O₂ failure:
N₂O supply is blocked and N₂O can no longer be
used.
- If there is an N₂O or air shortage or failure:
no alarm. The float in the right-hand flowmeter tube
falls, or drops right down to zero.

The cause of a disruption in supply can often be rectified quickly:

- Probe disconnected from medical pipeline system or
in standby position? Push probe in as far as it will go.
- Supply hose kinked? Straighten hose.
- Cylinder empty? Replace with full one.

If the cause cannot be rectified quickly:

- If there is an O₂ failure (medical air is not available):
ventilate patient manually with, for instance,
Resutator 2000.
- If there is an O₂ and/or N₂O failure:
use emergency supply from portable cylinder pack.

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*) Not allowed according to the European Standard.

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Additional instructions

- for Titus with S-ORC
- for Titus with air supply

Titus with S-ORC

S-ORC is a flow-dependent safety device.

Depending on the O₂ flow through the measuring tube, the maximum N₂O flow is limited so that the O₂ percentage in the fresh gas cannot fall below a minimum value (21 vol. %)

This means that N₂O supply is blocked not only when there is an O₂ failure, but also if the O₂ flow valve has been accidentally closed.

- Open O₂ flow valve.

Titus with air supply

If there is an O₂ failure this machine switches automatically to air supply, whatever the position of the change-over switch below the flow valves. The N₂O flow is blocked and air now flows through the N₂O flowmeter tube.

- Ventilate patient with Air.

If the patient was being ventilated with O₂ and N₂O anaesthetic gases before the O₂ failure occurred:

- use alternative or supplementary anaesthesia methods; maintain depth of anaesthesia.

If there is an air failure:

no alarm; float in the N₂O/Air flowmeter tube drops to zero. O₂ can still be administered.

After each O₂ failure:

Reset O₂ shortage alarm

Once O₂ supply is re-established, there must be a supply pressure of at least 2.7 bar at the shortage alarm for at least 20 seconds. During this period, do not switch on any appliances which use O₂, such as O₂ flush or O₂ fresh gas flow.

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Continuous fresh gas flow through the CO_2 absorber dries out the soda lime.

Does the moisture fall below minimum, undesired reactions, independent of the kind of lime and the used inhalation anaesthetic agent, arise:

- reduced CO_2 absorption,
- formation of CO,
- absorption and decomposition of the inhalation anaesthetic agent,
- increased heat development in the absorber and with this increased respiratory gas temperature.

These reactions can endanger the patient in form of CO poisoning, insufficient depth of anaesthesia and burns of the respiratory tracts.

N_2O is detrimental to the environment!

Do not use more fresh gas than necessary!

- Turn rotary knob on Ventillog to »0« position.
- Put O_2 , N_2O and air probes in »standby position« or disconnect.

When supply is from cylinders:

- Close cylinder valves.
- To relieve pressure on the anaesthetic machine, open N_2O flow valve first, then open O_2 flow valve until the floats in the flowmeter tubes have dropped completely.
- Close flow valves again.
- Switch off monitors.

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Care


Stripping down machine

Ancillary equipment – such as Circle System, Ventillog, monitors, vapor, etc. – see relevant Instructions for Use.

- Disconnect all probes from wall terminals – the hoses may remain on the machine.
- Disconnect fresh gas hose from Circle System – the hose remains on the anaesthetic machine.
- Pneumatic switching valve: loosen all hoses, disconnect switching valve from Circle System and also loosen hoses on the ventilator.
- Strip down Circle System, remove from anaesthetic machine and **proceed according to relevant Instructions for Use**.

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Magill Inhalation Device

- 1 Detach mask from mask nozzle, pull mask nozzle off the breathing hose, disconnect breathing hose and 2.3 L breathing bag from connection.
- 2 Unscrew connector from common-gas outlet.

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Care list for Titus anaesthetic apparatus

Only for non-infectious patients.

When used with infectious patients, all parts in contact with breathing gas must additionally be sterilized after disinfection and cleaning.

What	How often			
Components which can be conditioned	Conditioning intervals ¹⁾			
	with Filter on Y-piece	with Filter on inspiration and expiration port	without Filter	
Titus Basic				
Power cable, compressed gas hoses, equipotential bonding lead		daily	monthly	
Anaesthetic vapouriser				
Vapor 19.n / Vapor 2000		daily		
Circle system 9				
Metal parts	weekly	weekly	daily	
CO2-absorber with insert	weekly	weekly	daily	
ISO changeover valve	weekly	weekly	daily	
FGE changeover valve	weekly	weekly	daily	
Valve disks	weekly	weekly	daily	
Breathing bag with breathing hose connector	daily	daily	per patient	
Breathing hoses	daily	per patient	per patient	
Masks	per patient	per patient	per patient	
Fresh gas hose	weekly	weekly	weekly	
Y-Piece	daily	per patient	per patient	
Water trap	daily	per patient	per patient	
Pressure gauge hose	weekly	weekly	daily	
Probe scavenging hose or Probe return hose		weekly		
Monitor				
PM 8050, PM 8060	weekly	weekly	weekly	
Volumeter 3000	weekly	weekly	daily	
Volumeter 2000 K (for children)	weekly	weekly	daily	
Oxidig	daily	daily	daily	
Flow-sensor	weekly	weekly	daily	
O2-sensor	weekly	weekly	daily	
Temperature sensor	weekly	daily	daily	
Water trap (monitoring)	daily	daily	daily	
Filter				
Filter of the measured gas recirculation system	weekly	weekly	daily	
Microbial filter 645 St	-	daily	-	
AGS-System				
Scavenging hose with connector		monthly		
AGS housing		weekly		
Container for AGS buffer volume		weekly		
AGS transfer hose		weekly		
AGS flow tube (without filter)		weekly		
Waste gas connector with tube		weekly		
Emergency breathing				
Emergency breathing bag, its diaphragm and upper and lower valve parts (Dräger Resultator 2000)		after each use		
silicone sleeve of secretion collecting bottle and rinsing bottle, their lid with float, suction hose and inspection window		as requirement, but at least daily		

1) The conditioning intervals depend on the use and position of the filters. the table is merely intended as a rough guide. The instructions of the hospital's hygiene officer shall prevail!

2) Only use cleaning agents, no alkaline or chlorine-releasing substances. Risk of corrosion.

3) Surface disinfectants based on aldehydes, alcohols and quaternary ammonium compounds, such as Buraton 10 F (Schülke & Mayr GmbH) or Incidin Perfekt (Henkel Hygiene GmbH) can be used to ensure compatibility.

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4) Do not use any agents containing alcohol.

5) Do not submerge in water.

6) Disinfect the flow sensor in 70%ethanol solution for approx 1 hour and leave to dry in air for at least 30 minutes.

7) Use only pH-neutral detergent and distilled water.

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Cleaning and disinfection

- Wipe off heavy dirt with disposable cloths.

For disinfecting, use preparations of the surface disinfectants group – to avoid damaging the materials, use preparations based on the following active substances:

- Aldehydes
- Quaternary ammonium compounds.

Do not use the following products, because they may damage the material:

- Preparations based on halogen-separating compounds
- Strong organic acids
- Oxygen-separating compounds.

Recommendation for users in Germany:

- Disinfectants entered in the DGHM list – this list also specifies the active substance on which each disinfectant is based –
DGHM – Deutsche Gesellschaft für Hygiene und Mikrobiologie – DGHM-List available from mhp-Verlag GmbH, Wiesbaden.
For countries where the DGHM list is not known, the active bases specified above are recommended.

Wipe disinfecting

Titus, Vapor, Ventillog, medical gas connection hoses and control line of the changeover valve:

- Do not use any cleaning agents containing alcohol –
recommendation: e.g. Buraton 10 F –
Schülke & Mayr, Germany –
observe the manufacturer's Instructions for Use.

In cleaning and disinfecting machines

disinfect in liquid medium at high temperature

Connecting hose, bellows, patient system, FGE changeover valve, changeover valve and waste gas nozzle –

- 93 °C – 10 minutes –
- only with cleaning agent – do not add disinfectant – parts may corrode!
- **After disinfection in washing machine, always sterilise the FGE changeover valve at 134 °C in autoclave.**
Otherwise, operation may be disrupted because residual liquid may remain in the control channels!

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Sterilisation

Steam autoclaving (see appropriate column in Table, page 49).

Steam autoclaving accelerates the natural ageing process of rubber parts. Frequent checks for leaks and firm connection are essential.

Sterilise microbial filter in steam at 134 °C only.
Multiple sterilisation is possible (see relevant Instructions for Use).
Mark each sterilisation on the filter.

After care:

Assemble the device and check readiness for operation as described on page 52ff!

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Service/Maintenance intervals

Clean and disinfect device or device parts before each maintenance operation, including when returning for repair!

Apparatus Inspection	Every 6 months by an expert
Safety inspection	Every 6 months by an expert
O ₂ Sensor	Replace on reaching expiry date or when calibration is no longer possible
Control membrane of the O ₂ shortage threshold signal	Should be replaced every 4 years by an expert
Membrane of S-ORC	Should be replaced every 5 years by an expert
Membrane of the pressure reducer	Should be replaced every 6 years by an expert

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Testing Readiness for Operation

– after completion of each care procedure

Completely re-assemble the machine and supply with compressed gases (cylinders – if available – filled); CO₂ absorber filled.

O₂ cylinders are completely charged if the O₂ pressure gauge indicates 200 bar. At this pressure, a 3-L cylinder contains 600 litres and 11-L cylinder contains 2200 litres of expanded gas.

The N₂O pressure amounts to about 50 bar, as long as liquid nitrous oxide is available in the cylinder. When gas is extracted, the cylinder pressure may temporarily drop below 50 bar, as a function of the related temperature reduction. The liquid content of nitrous oxide cylinders can only be determined by weighing. A charged 3-L cylinder (2.25 kg charging weight) contains about 1125 litres, and a charged 11-L cylinder (8 kg charging weight) contains about 4000 litres of expanded nitrous oxide.

Checking compressed gas connections

Check all supply hoses (from pipeline system or cylinders) for firm connection, leaks and general condition (visual check).

Checking safety features

(Description of safety features, page 61).

For all devices – with or without S-ORC, and with or without an additional air supply:

- N₂O must not flow if the O₂ compressed gas supply has been interrupted or has failed
 - = effect of N₂O cut-out for machines without S-ORC
 - = effect of S-ORC if machine is equipped with it (see label »S-ORC« next to the N₂O measuring tube).
- If there is an O₂ shortage or failure, an acoustic alarm sounds.

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Checking O₂ shortage alarm and N₂O cut-out (all Titus models)

- For Titus Air: set switch on flowmeter block to »N₂O« position.
- Open O₂ and N₂O; flow approx. 3 L/min. ● Cut O₂ supply – disconnect probe from wall terminal or close cylinder valve; after a short time an alarm tone must sound for at least 7 seconds and the N₂O flow must return to zero.
- Close both flow-control valves and re-establish O₂ supply. After 20 seconds at a minimum of 2.7 bar operating pressure, the O₂ shortage alarm will be operational again. During this time, no O₂ or N₂O should be delivered by the machine.

Checking S-ORC (machines with S-ORC)

S-ORC limits the max. N₂O flow as a function of the O₂ flow, and thereby guarantees a minimum O₂ concentration. S-ORC also prevents N₂O flow if the O₂ valve is closed.

- Set O₂ flow to 1 L/min, and then open the N₂O flow control valve completely. There should be an O₂ concentration of between 21 vol.% and 27 vol.%, or an N₂O flow of between 2.7 L/min and 3.7 L/min.
- Close the O₂ flow control valve: the N₂O flow must return to zero, i.e. there must be an O₂ concentration of 100 vol.%.
- Open the O₂ valve again; N₂O flows again.
- Check whether the flow can be set to 10 L/min for O₂, N₂O and Air.

Checking selector switch and gas type


only on machines which have an air supply.

Machines without air supply; carry out the check with 100 vol.% O₂ and with O₂/N₂O mixture 1:1.

- Keep O₂ measuring device (e.g. Oxydig) with sensor in operation for the test.
- Unscrew hose from fresh-gas outlet on Titus.

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**Selector switch in »N₂O« position**

- Set O₂ flow = 3 L/min.
- Hold O₂ sensor under the fresh gas outlet and allow gas to flow over the sensor: O₂ concentration on the O₂ measuring device must clearly increase towards 100 vol.% O₂.
- Open N₂O flow valve, set flow = 3 L/min: the O₂ concentration must clearly drop towards 50 vol.%.

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Selector switch to »Air« position

- Open both flow control valves and set the same flow – e.g. 3 L/min – on both flowmeter tubes.
 - Test with O₂ measuring device: The display should read about 60 vol.% O₂.
 - Interrupt the O₂ supply. The O₂ shortage alarm must sound and the air-flow must remain unchanged.
 - Test with O₂ measuring device: The display should read 21 ±3 vol.% O₂.
 - Close both flow valves; reconnect O₂ supply.
- Keep fresh-gas outlet open for the next test:

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
Finally:

- Reconnect the fresh gas hose to the common-gas outlet on the anaesthetic machine. Replace Vapor.

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O₂ Flush (Bypass)

- Press O₂-Flush: the Y-piece on the circle system must emit a strong flow of gas – do not close Y-piece. Release button. The O₂ flow should cease.

Circle System 9

- Check that the fresh gas hose and its connections are in perfect condition (visual check).
- Check functioning according to »Circle System 9« Instructions for Use.

Ventilog 2 or 3

- Ventilog supplied with compressed gas, hose in order.
- Check all connections between Ventilog, pneumatic switching valve, circle system and breathing bag for firm connection.
- Check functioning of Ventilog according to relevant Instructions for Use.

Bronchial Aspirator

Check functioning according to relevant Instructions for Use.

Microbial Filter

Check condition and fitting of filter in accordance with the recommendations contained in the »Microbial Filter« Instructions for Use.

Vapor

- Only plug-in Vapor: Check that the Vapor plug-in adaptor is fitted level onto the mounting (on the anaesthetic machine), and that the system is locked.
- Check filling level; refill if necessary.
- Check functioning of Vapor according to its Instructions for Use.

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Anaesthetic gas scavenging system (AGS)

- Check functioning according to the relevant Instructions for Use.

Anaesthetic gas scavenging*

For scavenging via ejector unit:

- Check that the circle system and ventilator each have a connecting hose (see page 22).
- Connect hoses to T-piece.
- Connect probe on anaesthetic scavenging hose into terminal unit on the scavenging unit – »green« indicator on the scavenging terminal unit shows that the unit is switched on.

If no ejector unit is available:

- Connect an anaesthetic filtering device* to each of the waste gas sockets – replace the filters!
- Check that the filters are firmly fitted to in rubber cuffs.

Manual ventilation bag for emergency ventilation

- Check functioning by pumping by hand: when the bag is squeezed, the air escaping from the mask cone must be audible and obvious. When the bag is released, it must rapidly resume its former shape.
- When mask connector is closed (with palm of hand), it should only be possible to squeeze the bag a little.

Monitors

- Check functioning of monitors according to relevant Instructions for Use.

After these tests have been carried out successfully, the machine is ready for operation.

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Functional Diagram

The diagram on page 60 applies to Titus with its maximum configuration, i.e. **with S-ORC** (In the area of validity of the European Standard only devices with S-ORC are licensed.) and **with Air** as a third gas.

The machine is supplied with O₂, N₂O and air from a pipeline system (PS) via wall terminal units. Alternatively, O₂ and N₂O may be supplied from compressed gas cylinders **18**. Cylinder pressure is displayed on a pressure gauge **17** and reduced by a pressure reducer to 5 bar.

Pressure reducers **5.1** reduce the O₂ and N₂O pressure to 2 or 3 bar; pressure reducer **5.2** reduces the air pressure from 5 to 1.5 bar.

N₂O and air can now be administered alternately, depending on whether the selector switch **6** is in the N₂O or »Air« position. Both gases are administered via the same (right) flowmeter tube (see also the diagram in the Technical Data).


The O₂ pressure is monitored by an O₂ shortage alarm **2, 3, 4** (with whistle): the alarm sounds when the O₂ pressure is lower than 2.2 bar. If O₂ pressure continues to drop, the N₂O supply – whatever the position of the selector switch **6** – is interrupted and the machine switches automatically to »Air« via the control block **7**. If the selector switch **6** is in the »Air« position, medical air remains switched on even during O₂ shortage. If the O₂ flow is inadequate, the S-ORC **8** limits the N₂O flow, so that O₂ concentration in the mixed gas does not drop below 21 vol.%. If the O₂ flow is switched off, or the O₂ supply interrupted, the N₂O flow is stopped.

O₂, N₂O and air flows can be set on the flow control valves **9** and readings taken on the flowmeters **10**.


Downstream of the flowmeters, the gases converge and flow through connections **12** or **13** to the anaesthetic vaporiser **14**, where the anaesthetic agent is added. When a detachable plug-in vaporiser **14** is being used, the probe system **13** blocks off the connections to the outside when the vaporiser is detached and opens an internal connection so that the gas can flow to the fresh-gas outlet **16**.

The safety valve **11** (only on machines with S-ORC), upstream of the inlet to Vapor, limits the pressure in the flowmeter to about 0.9 bar.

With the O₂ flush **15** (bypass), an O₂ flow of about 55 L/min (depending on supply pressure) can be added to the fresh gas, without affecting the pressure ratios.

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- | | | | |
|-----|---|----|--|
| 1 | Compressed gas connections with filter | 10 | Flowmeters |
| 2 | Non-return valve | 11 | Safety valve (only for S-ORC) |
| 3 | Reservoir | 12 | Vapor, fixed |
| 4 | Alarm for O ₂ shortage | 13 | Plug-in system |
| 5.1 | Pressure reducers 5/3 bar for O ₂ and 5/2 bar for N ₂ O;
5/1.5 bar on machines without S-ORC | 14 | Vapor, plug-in type |
| 5.2 | Pressure reducer 5/1.5 bar (for Air) | 15 | O ₂ flush (bypass) |
| 6 | Selector switch »N ₂ O«/»Air« | 16 | Fresh-gas outlet |
| 7 | Control block | 17 | Pressure gauge for cylinder pressure with 200/5 bar pressure
reducer. |
| 8.1 | S-ORC control element for minimum O ₂ concentration | 18 | O ₂ and N ₂ O compressed gas cylinders. |
| 8.2 | Flow resistance | | |
| 9 | Flow valves | | |

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All Titus models:

O₂ shortage threshold

Alarm limit	Warning at the latest when O ₂ pressure has fallen below 1.8 bar.
Warning	At least 7 seconds of continuous sound; alarm tone cannot be switched off.
O ₂ pressure necessary to reset the alarm	≥2,7 bar for at least 20 seconds; do not switch on any appliances which use O ₂ during this period such as O ₂ flush or O ₂ flow.

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
Titus with S-ORC:

Oxygen Ratio Controller

To ensure that the O₂ concentration in the fresh gas always remains above 21 vol.%, the anaesthetic machine is equipped with an O₂ ratio controller. S-ORC limits the N₂O flow as a function of the O₂ flow. Up to an O₂ flow of 200 mL/min, the N₂O flow is stopped, and above 200 mL/min it is slowly released. For an O₂ flow above 300 mL/min, N₂O can be added while the S-ORC prevents the O₂ concentration falling below 21 vol.%.

The graph shows the range of adjustable O₂/N₂O ratios.

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S-ORC is not oxygen monitoring

By adding anaesthetic agent, the O₂ concentration in the fresh gas is reduced.

Therefore:

Check the O₂ concentration with monitor!

For Titus with »Air« supplied as the third gas control of fresh gas flow – one of the gases fails – depends on the position of the selector switch (»N₂O« or »Air«).

Switch in »N₂O« position	O ₂ + N ₂ O can be metered; air supply is automatically cut off.
If O ₂ is inadequate or fails	Automatic switch-over to air; N ₂ O is cut off. The selector switch remains unchanged in position in the »N ₂ O« position.
If N ₂ O fails	O ₂ can still be metered. No alarm; float in N ₂ O flowmeter tube drops to zero.
Switch in »Air« position	O ₂ + N ₂ O can be metered; air supply is automatically cut off.
If O ₂ is inadequate or fails	S-ORC is not effective in the »Air« position. Only air can be metered without limit.
If Air fails	O ₂ can still be metered. No alarm; float in N ₂ O flowmeter tube drops to zero.

Titus without S-ORC and without air supply

If O ₂ is inadequate or fails	Acoustic alarm, the N ₂ O cut-out device switches off the N ₂ O supply.
If N ₂ O fails	O ₂ can still be metered. No alarm; float in N ₂ O flowmeter tube drops to zero.

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Ambient Conditions

During operation

Temperature	15 to 35 °C
(machine temperature must be at ambient temperature)	
Air pressure	700 to 1060 hPa
Humidity	20 to 80 % rel. humidity

During storage

Temperature	–10 to 60 °C
Air pressure	500 to 1060 hPa
Humidity	0 to 80 % rel. humidity

Noise emission

During normal operation	none
Alarm tone	approx. 73 dB (A)

Gas supply

From pipeline system (PS) with the following requirements:

O ₂ pressure	2.7 to 6.0 bar
N ₂ O pressure	2.7 to 6.0 bar
Air pressure	2.7 to 6.0 bar
Gas purity	in accordance with requirements in European Pharmacopeia.

Connectors on machine (conform to EN 739)

for O ₂	NIST
for N ₂ O	NIST
for medical air	NIST

Optional O₂ and N₂O supply from
3-L or 11-L cylinders

Gas purity:	in accordance with requirements in European Pharmacopeia.
Pressure reduction:	via pressure reducer on cylinders, outlet pressure 5 bar.

Gas dosage

● O ₂ dosage	Fine dosing valve (5 revolutions)
Measuring range	0,5 to 10 L/min O ₂ (one flowmeter tube)
Tolerance (20 °C, 1013 hPa)	±10 % of displayed value +15, –5 % for lowest scale value

Flowmeters for low-flow (O₂):

1st flowmeter tube	
Measuring range	(0.02 to 0.5) L/min

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Flowmeters for low-flow (N₂O):

1st flowmeter tube	
Measuring range	(0.02 to 0.5) L/min
Tolerance	±10 % or for: lower scale values -20 % / +10 % upper scale value ±5 % of the displayed value,
2nd flowmeter tube	
Measuring range	(0.55 to 10) L/min
Tolerance	±10 % or for: lower scale values +20 % / -5 % of the displayed value,
If Titus is supplied with additional air:	
● »Air«/»N ₂ O« selector switch, manually operated	
»Air« position:	O ₂ + Air can be dosage.
»N ₂ O« position:	O ₂ + N ₂ O can be dosage.
● Air metering	with N ₂ O dosing valve, Air flows through the N ₂ O flowmeters.
Measuring range	0.55 to 10 L/min Air *) Display as shown in graph
Tolerance (20 °C, 1013 hPa)	±20 % of the displayed value

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*) In low-flow version (double flowmeter tubes)

1st flowmeter 0.02 to 0.5 L/min
2nd flowmeter 0.55 to 10 L/min

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
O₂ Flush (Bypass)

Supply volume at 5 bar: approx. 55 L/min O₂
Supply volume at 2.7 bar: approx. 35 L/min O₂

Anaesthetic agent dosage with:

- Dräger-Vapor® 19.n; Setting range for
Halothane-Vapor: 0.2 to 4 vol.% Halothane
Enflurane-Vapor: 0.2 to 5 vol.% Enflurane
Isoflurane-Vapor: 0.2 to 5 vol.% Isoflurane
Sevoflurane-Vapor: 0.3 to 8 vol.% Sevoflurane
- Dräger-Vapor® 2000; Setting range for
Halothane-Vapor: 0.2 to 6 vol.% Halothane
Enflurane-Vapor: 0.2 to 8 vol.% Enflurane
Isoflurane-Vapor: 0.2 to 6 vol.% Isoflurane
Sevoflurane-Vapor: 0.3 to 8 vol.% Sevoflurane
- Titus with plug-in Vapor:
Rapid change system for Vapor.
The connections are automatically made and sealed
against atmosphere when Vapor is removed.
- Equipment and technical data of Vapor:
see "Vapor 19.n or Vapor 2000
Instructions for Use"

Dosage of air with the N₂ O flowmeter tube

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Common-gas outlet on Titus

- Connector according to ISO:
- Straight nozzle, outlet pointing horizontally to left

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- O₂-Mangelsignal
- N₂O-Sperre (wenn ohne S-ORC)
- S-ORC (Sensitive Oxygen Ratio Controller), N₂O-Sperre entfällt.

Safety features

Description see page 61

- O₂ shortage alarm
- N₂O cut-out device (if no S-ORC)
- S-ORC (Sensitive Oxygen Ratio Controller), with S-ORC the N₂O cut-out device is not needed.

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Circle system

Equipment and technical data: see »Circle System 9« Instructions for Use

Fresh gas inlet connector:
male thread M 16 x 1.5 mm

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Dimensions and weights

Weight:


Machine with Vapor 19.n
and circle system 22 kg*

Machine with Vapor 19.n
without circle system 16 kg*

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*) The weight increases by around 1 kg at the use of the Vapor 2000.

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These Instructions for Use apply only to
Titus
with Serial No.:

If no Serial No. has been filled in by Dräger
these Instructions for Use are provided for
general information only and are not
intended for use with any specific
machine or device.

Directive 93/42/EEC
concerning Medical Devices

Dräger Medical AG & Co. KGaA

🏠 Moislinger Allee 53 – 55
D-23542 Lübeck
☎ +49 451 8 82 - 0
FAX +49 451 8 82-20 80
💻 <http://www.draeger.com>

Dräger Medical AG & Co. KGaA

Germany
🏠 Moislinger Allee 53 – 55
D-23542 Lübeck
☎ +49 451 8 82 - 0
FAX +49 451 8 82-20 80
💻 <http://www.draeger.com>